BabyTemp Schnuller-Fieberthermometer

DE

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Art.Nr.: 9633

Model: TM03

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Model: TM03

Special features:

The pacifier thermometer is a maximum thermometer, i.e. the highest temperature is displayed after only 90 seconds. The most recent measured value is always saved.

Safety instructions:

- The thermometer is not suitable for use as a pacifier.
- The thermometer is suitable exclusively for oral measurement of the body temperature of babies.
- The thermometer can not be opened and the battery cannot be changed for safety reasons.

Operating manual:

Press the I/O button to switch on. A short signal tone shows that the device is in the „switched on“ condition. A segment test takes place simultaneously (checks all display segments and carries out an internal function test (see under „Function test“)).

After this, if the temperature measured by the temperature sensor is less than 32°C, the display shows an „Lo“ and „°C“. After around one second, the „°C“ starts to blink. The „Lo“ stands for low (less than 32°C).

The current temperature is continuously displayed during measurement, and the „°C“ blinks during the entire measurement time. The measurement ends as soon as temperature stability is achieved. An acoustic signal is produced (10 x 3 quick beeps).

After the measurement is completed, the temperature value measured remains in the display until the pacifier thermometer is switched off. To turn the device off after use, press the I/O button briefly. The signal tone will again be sounded. If the device is not switched off, a signal tone is sounded after around 9 minutes and the device then switches off automatically.

Display of last measured value:

If you press the I/O , the last measured value will be shown in the display.

Function test:

The thermometer is checked for correct functioning automatically every time it is switched on.

After display of last measured value, the screen shows 37,0 °C if correct operation.

Important:

The readings displayed by the thermometer may be too low under some conditions (e.g. if the child is crying or persistently breathing through his mouth). Therefore, in case of doubt, or if there are signs of fever (reddening of skin, sweating, agitation, unwillingness to drink), please always perform a follow-up rectal measurement if the result is more than 38,5°C.

Cleaning and disinfection:

Please only use an alcohol pad. Each pharmacy sells alcohol pads. We recommend the use of Mucocit®-F2000 or Helpipur®Hplus to disinfect the device during use in medical facilities. Both are effective against all hospital-typical pathogens and can be purchased from any dispensing chemist.

Battery symbol:

If the battery symbol  displayed, the life time of the thermometer is over. The battery cannot be changed for safety reasons.

Measurement checks:

Checking the accuracy of the pacifier thermometer is not necessary in most cases. If you require a measurement check (at your expense), please send the pacifier thermometer to our customer service address. Please make sure it is well-packaged and that the correct postage stamps are applied.

Guarantee:

We give a guarantee period of 2 years from the date of purchase. We guarantee perfect design of this article. We will rectify all material or manufacturing errors during the guarantee period. If, contrary to expectations, defects are found, please send the device back to our customer service department. Please enclose the proof of purchase and a description of the defect so that the problem can be rectified as soon as possible.

The guarantee excludes damages caused by improper handling, related to the battery or to the packaging.

Disposal note:

The disposal of batteries is regulated. Dispose of the thermometer as prescribed, i.e., at a collection point for electrical scrap.

Specifications:

Measuring range:	32.0°C to 42.9°C Temperature below 32.0°C: display "Lo" for low. Temperature above 43.0°C: display "Hi" for high.
Measurement accuracy:	0.1°C between 35.5°C and 42.0 °C at a surrounding temperature of 20°C to 24°C.
Self-test:	Automatic internal check of test value 37.0°C.
Display:	3-digit liquid-crystal (LCD display).
Battery display:	The  shows when the battery is discharged.
Storage temperature:	0°C to 50°C.
Changes in terms of technical progress reserved.	

Legal basis for marketing the device in the EU:

The manufacturer was certified according to the guidelines of 93/42/EEC for medical products. This device complies with the EN60601-1-2 European standard on electromagnetic compatibility and meets the requirements of EN 12470-3.

The CE marking certifies conformity with the basic requirements of the 93/42/EEC guidelines for medical products.

Explanation of symbols:

	Manufacturer	 Cotronic Technology Ltd Floor 4~6,Block 7, west of ZhouShi Road, XiXiang Street, Bao An Zone, ShenZhen, GuangDong, China.
	Read the instruction manual	
	Product type BF	
	Disposal notes	
	Production period (YYYY/MM)	
	Die CE-Kennzeichnung der Benannten Stelle	EC REP Company Address WellKang Ltd., Enterprise Hub, NW Business Complex, 1 Beraghmore Road,Derry, BT48 8SE, Northern Ireland, UK

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reer GmbH • Mühlstr. 41 • 71229 Leonberg • Germany • www.reer.de