BOUTEILLE DE RECUEIL DE LAIT MATERNEL 180ML (X1)
BREAST MILK COLLECTION BOTTLE 180ML (X1)

- FR** Notice d'utilisation
- EN** Instructions
- NL** Handleiding
- DE** Gebrauchsanweisung
- ES** Folleto de Instrucciones
- IT** Istruzioni per l'uso
- PT** Instruções de uso
- PL** Instrukcja użycia
- RU** Инструкция
- CN** 操作说明
- CT** 使用說明
- AR** تعليمات الاستخدام

FR

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EN

Read carefully and keep for future reference.

This bottle is designed to collect milk obtained using your Béaba hands-free electric breast pump. You can also store the milk in it by using the cap provided. The two elements (bottle and lid) can be used in the fridge and freezer.

WARNING!

- Always remove the bottle from the motor unit before storing.
- Always ensure that the bottle and the lid are not cracked or damaged before each use.
- The bottle must not be used to reheat the breast milk. Do not use a microwave or bain-marie to reheat milk.
- Do not defrost the milk in the open but in the fridge.
- Wash your hands carefully before touching the bottle and avoid touching the inside.
- Never refreeze defrosted breast milk.
- Always check the temperature of the milk before serving it to baby.

1. Preparation

- **Cleaning:** Clean all the parts of the bottle (bottle, lid) carefully in warm soapy water before and after each use and rinse well. Leave to dry in the open and check that they are dry before using.
- **Before using for the first time:** Sterilisation: After washing them, sterilise the bottle and lid by boiling them for 5 minutes or using a steam steriliser.

2. Instructions for use:

Wash your hands carefully before commencing.

- **Collecting milk:**
 - Connect the bottle to the breast pump following the instructions provided with your Béaba breast pump.
 - Make sure that the connection is correctly sealed to avoid any leaks.
 - Start expressing your milk following the instructions provided with your Béaba breast pump.

- **If you choose to store your milk in the bottle:**

Remove the bottle from the electric unit after expressing your milk and close it with the lid supplied.

Remember to label the bottle with the date and time the milk was expressed and the volume of milk collected.

Be sure not to fill the bottle to the maximum limit (freezing will increase the volume of milk in the bottle).

Breast milk can be refrigerated, frozen and defrosted in the BEABA bottle. Please follow the instructions below:

Breast milk	Storage	Shelf life
Fresh	Room temperature (22°C max)	4 to 6 hours
	Fridge (0 to 4°C)	3 to 5 days
	Freezer (freezer compartment of fridge)	2 weeks
	Freezer (-18°C)	6 to 9 months
Defrosted	Room temperature (22°C max)	2 hours
	Fridge (0 to 4°C)	24 hours

- **Do not refreeze after defrosting.**
- **Please note that these times cannot be combined.**
- These storage and defrosting rules for breast milk are recommendations only.
- **Consult a healthcare professional or a breastfeeding specialist for further information.**

• **Defrosting milk:**

- Do not use a microwave or bain-marie to defrost breast milk.
- Defrost the milk gradually in the fridge.
- In the event of an emergency, you can defrost breast milk by placing the bottle in cold or lukewarm water (max. 37°C) and using it immediately.
- To reheat the milk, pour it into a suitable container and refer to the manufacturer's instructions.

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BEABA ASIA LTD